



# Recombinant Mouse Tumor necrosis factor receptor superfamily member 11A (Tnfrsf11a)

<b>Product Code</b>	CSB-CF023968MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O35305
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	VTPPCTQERHYEHLGRCCSRCEPGKYLSSKCTPTSDSVCLPCGPDEYLDWTN EEDKCLLHKVCDAGKALVAVDPGNHTAPRRCACTAGYHWNSDCECCRRNTE CAPGFGAQHPLQLNKDTVCTPCLLGGFFSDVFSSTDKCKPWTNCTLLGKLEAH QGTTESDVVCSSTMLRRPPKEAQAYLPSLIVLLLIFISVVVVAIIIFGVYYRKGK KALTANLWNWVNDACSSLSGNKESGDRCAGSHSATSSQQEVCEGILLMTR EEKMVPEDGAGVCGPVCAAGGPWAEVRDSRTFTLVSEVETQGDLSRKIPTED EYTD RPSQPSTGSLLLIQGSKSIPPFQE PLEVGENDSL SQCFTGTTESTVDSE GCDFT EPPSRTDSMPVSPEKHLTKEIEGDSCLPWWVSSNSTDGYTGSGNTPG EDHEPFGSLKCGPLPQCAYSMGFSEAAASMAEAGVRPQDRADERGASGS GSSPSDQPPASGNVTGNSNSTFISSGQVMNFKGDIIVVYVSQTSQEGPGSAE PESEPVGRPVQEETLAHRDSFAGTAPRFPDVCATGAGLQEQQGAPRQKDGTS RPVQEQQGAQTSLHTQGSGQCAE
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Tnfrsf11a
<b>Protein Names</b>	Recommended name: Tumor necrosis factor receptor superfamily member 11A Alternative name(s): Osteoclast differentiation factor receptor Short name= ODFR Receptor activator of NF-KB CD_antigen= CD265
<b>Expression Region</b>	31-625
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.