



Recombinant Mouse Interleukin-27 receptor subunit alpha (Il27ra)

Product Code	CSB-CF011644MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70394
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	TRPHGSPGPLQCYSVGPLGILNCSWEPLGDLETTPVLYHQSQKYHPNRVWEV KVPSKQSWVTIPREQFTMADKLLIWGTQKGRPLWSSVSVNLETQMKPDTQPQIF SQVDISEEATLEATVQWAPPVWPPQKVLICQFRYKECQAETWTRLEPQLKTD GLTPVEMQNLEPGTCYQVSGRCQVENGYPWGEWSSPLSFQTPFLDPEDVW VSGTVCETSGKRAALLVWKDPRPCVQVTTYVWFGAGDITTTQEEVPCCKSPV PAWMEWAVVSPGNSTSWVPPTNLSLVCLAPESAPCDVGVSSADGSPGIKVT WKQGTRKPLEYVVDWAQDGDSDLKLNWTRLPPGNLSTLLPGEFKGGVPYRIT VTAVYSGGLAAAPSVWGFREELVPLAGPAVWRLPDDPPGTPVVAWGEVPRH QLRGQATHYTFICIQSRGLSTVCRNVSSQTQTATLPNLHLGSFKLWTVSTVAG QGPPGNLSLHLPDNRIRWKALPWFLSLWGLLLMGCGLSLASTRCLQARCLH WRHKLLPQWIWERVPDPANSNSGQPYIKEVSLPQPPKDGPILEVEEVELQPVV ESPKASAPIYSGYEKHFLPTPEELGLLV
Source	in vitro E.coli expression system
Target Names	Il27ra
Protein Names	Recommended name: Interleukin-27 receptor subunit alpha Short name= IL-27 receptor subunit alpha Short name= IL-27R subunit alpha Short name= IL-27R-alpha Short name= IL-27RA Alternative name(s): Type I T-cell cytokine receptor Short nam
Expression Region	25-623
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.