



Recombinant Human Leukocyte immunoglobulin-like receptor subfamily B member 3 (LILRB3)

Product Code	CSB-CF012941HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75022
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	GPFPKPTLWAEPGSVISWVSPVTIWCQGSQEAQEYRLHKEGSPEPLDRNNPL EPKNKARFSIPSMTEHHAGRYRCHYYSSAGWSEPSDPLEMVMGTGAYSKPTLS ALPSPVVASGGNMTLRCGSQKGYHHFVLMKEGEHQLPRTLDSQQLHSRGFQ ALFPVGPVTPSHRWRFTCYYYTNTPWVWSHPSDPLEILPSGVSRKPSLLTLQ GPVLAPGQSLTLQCGSDVGYNRFVLYKEGERDFLQRPQGQPQAGLSQANFT LGPVSPSNGGQYRCYGAHNLSSEWSAPSDPLNILMAGQIYDVTSLSAQPGPT VASENVTLCCQSWWQFDTFLLTKEGAAHPPLRLRSMYGAHKYQAEFPMSP V TSAHAGTYRCYGSYSSNPHELLSHPSEPLELVVSGHSGGSSLPPTGPPSTPG LGRYLEVLIGVSVAFVLLL FLLL FLLL RRRQRHSHKRTSDQRKTD FQR PAGA AET EPKDRG LLRRSSPAADVQEENLYAAVKDTQSEDRVELDSQSPHDEDPQAVTY APVKHSSPRREMASPPSSLSGFLDTKDRQVEEDRQMDTEAAASEASQDVTY AQLHSLTLRRKATEPPPSQEGEPPAEPsiyatLAIH
Source	in vitro E.coli expression system
Target Names	LILRB3
Protein Names	Recommended name: Leukocyte immunoglobulin-like receptor subfamily B member 3 Short name= LIR-3 Short name= Leukocyte immunoglobulin-like receptor 3 Alternative name(s): CD85 antigen-like family member A Immunoglobulin-like transcript 5 Shor
Expression Region	24-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.