



Recombinant Acetolactate synthase-like protein (T26C12.1)

Product Code	CSB-CF011687CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O61856
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MVLFLLIAAIIIGLLLWKWLDVRSTDELTSMVKLLGSGNGQHVLSNAFQVDEKSK RHGGELVASVLKAHDVEEIFVLCGGHISPILVAAEKLGIKIVDTRHEVTAVFAAD AVARLRQSIGVAAVTAGPGLTNTITAVKNAQMAESPLLLIGGAAPTLLKGRGAL QDIDQMVLFRPLCKYVARVERLRDIVPTVREAIKAAKSGCPGPVFVEFPVDVLY PYELVVKEIGFNPNAKGFIQRALNFYL RCHVSRQFGNAWAPQTITPLPTNIPMP KSEKIQEIVQLVKSARKPVLLIGSQATLPPVKPADLVKAVEALGCPVFLGGMAR GLLGKDHPLQMRQVRRDALKDADLTILAGTVCDFRLSYGRTLSKKSIVALNR NSSQLTKNEKAFWNSDVSQADVATSLVQVANALGANHTTTPTEWVKSLREK DDEKESANAKKMEQKLTNGFLNPLNFLRTL DQSLPDDAILVADGGDFVGSAA YIVRPRGPLQWLDPGAFTLVGGGFALGAKTVYPKRPVYIIWGDGSCGYSLM EYDTFARHKLPVIGIVGNDACWTQIAREQVPMFQSSVAVDLARTRYDNVAKSL GSWGETIDESNADSARKILDEALAVCRSGEQSALVNVLIGKTDFFREGSISV
Source	in vitro E.coli expression system
Target Names	T26C12.1
Protein Names	Recommended name: Acetolactate synthase-like protein EC= 2.2.1.- Alternative name(s): IlvB-like protein
Expression Region	1-640
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.