



# Recombinant Human Mannosyl-oligosaccharide 1,2-alpha-mannosidase IA (MAN1A1)

<b>Product Code</b>	CSB-CF013383HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P33908
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	MPVGGLLPLFSSPAGGVLGGGLGGGGGRKGS GPAALRLTEKFVLLLVSFAFIT LCFGAIFFLPDSSKLLSGVLFHSSPALQPAADHKPGPGARAEDAAEGRARRRE EGAPGDPEAALEDNLRARENHERALREAKETLQKLPEEIQRDILLEKKKVAQD QLRDKAPFRGLPPVDFVPPIGVESREPADA AIREKRAKIKEMMKHAWNNYKGY AWGLNELKPI SKGGHSSSLFGNIKGATIVDALDTLFIMEMKHEFEEAKSWVEEN LDFNVNAEISVFEVNIRFVGGLLSAYYLSGEEIFRKKAVELGVKLLPAFHTPSGI PWALLNMKSGIGRNWPWASGGSSILAEFGTLHLEFMHLSHLSGNPIFAEKVM NIRTVLNKLEKPGQLYPNYLNPSSGQWGQHHSVVGGLGDSFYEYLLKAWLMS DKTDLEAKMYFDAVQAIETHLIRKSSSGLTYIAEWKGGLEHKMGHLTCFAG GMFALGADA APEGMAQHYLELGAEIARTCHESYNRTFMKLGPEAFRFDGGVE AIATRQNEKYYILRPEVMETMYMWRLTHDPKYRKWAVEAVEALENHCRVN GGYSGLRDVYLLHESYDDVQQSFFLAETLKYLYLIFSDDDLLPLEHWIFNSEAH LLPILPKDKKEVEIREE
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MAN1A1
<b>Protein Names</b>	Recommended name: Mannosyl-oligosaccharide 1,2-alpha-mannosidase IA EC= 3.2.1.113 Alternative name(s): Man(9)-alpha-mannosidase Short name= Man9-mannosidase Mannosidase alpha class 1A member 1 Processing alpha-1,2-mannosidase IA Short name= Alpha
<b>Expression Region</b>	1-653
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.