



Recombinant Rat Delta-like protein 1 (DII1)

Product Code	CSB-CF006947RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97677
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	QVWSSGVFELKLQEFVNKKGLLGNRNCCRGGSGPPCACRTFFRVCLKHYQA SVSPEPPCTYGSAVTAVLGVDSFSLPDGAGIDPAFSNPIRFPFGFTWPGTFLII EALHTDSPDDLATENPERLISRLTTQRHLTVGEEWSQDLHSSGRTDLRYSYRF VCDEHYGEGCSVFCRPRDDAFGHFTCGERGEKMCDPGWKGQYCTDPICLP GCDDQHGYCDKPGECKCRVGVWQGRYCDCEIRYPGCLHGTCQQPWQCNCQ EGWGGLFCNQDLNYCTHHKPCRNGATCTNTGQGSYTCSCRPGYTGANCEL EVDECAPSPCRNGGSCDLEDSYSCTCPPGFYGVKVELSAMTCADGPCFNG GRCSDNPDGGYTCHCPAGFSGFNCEKKIDLCSSSPCSNGAKCVDLGNSYLC RCQTGFSGRYCEDNVDDCASSPCANGGTCRDSVNDFSCTCPPGYTGRNCS APVSRCEHAPCHNGATCHQRGQRYMCECAQGYGGANCQFLLPEPPDLIVA AQGGSFPWVAVCAGVVLVLLLLLGC AAVVVCVRLKLQKHQPPDP CGGETET MNNLANCQREKDVSVSIIGATQIKNTNKKADFDHGDHGADKSSFKARYPTVDYN LIRDLKGDEATVRDAHSKRDTKCQSQGSVGEKSTSTLRGGVEVPDRKRPEVS YSTSKDTKYQSVYVLSAEKDECVIATEV
Source	in vitro E.coli expression system
Target Names	DII1
Protein Names	Recommended name: Delta-like protein 1 Alternative name(s): Drosophila Delta homolog 1 Short name= Delta1
Expression Region	18-714
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.