



# Recombinant Mouse Cyclic AMP-dependent transcription factor ATF-6 beta (Atf6b)

<b>Product Code</b>	CSB-CF002278MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O35451
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	AELMLLSEIADPTRFFFTDNLLSPEDWDSTLYSGLDEVAEEQAQLFRCVEQDVP FDSSSLDVGMDVSPPEPPWDPLPIFPDLQVKSEPSSPCSSSSLSSESSLSTE PPSQVPGVGEVLHVKMESLAPPLCLLGDDPASPFPETVQITVGSASDDLSDIQT KLEPASPSSSVHSEASLLSADSPSQPFIGEEVLEVKTESPSPPGCLLWDVPAS SLGAVQISMGSPDSSSGKAPATRKPLQPKPVVLTTPVPPRAGPTSAAVLL QPLVQQPAVSPVLIQGAIRVQPEGPAPAAPRPERKSIVPAPMPGNSCPPEVD AKLLKRQQRMIKNRESACQSRRKKKEYLQGLEARLQAVLADNQQLRRENAAL RRRLEALLAENSGLKLGSNGNRKVVCMVFLLFIAFNFGPVSISEPPPAPMSPRM SREEPRPQRHLLGFSEPGPAHGMEPLREAAQSPGEQQPSSAGRPSFRNLTA FPGGAKELLLRDLDQLFLSSDCRHFNRTESLRLADELSGWVQRHQRGRRKIP HRAQERQKSQRKKSPVVKPVPTQPPGPPERDPVGGQLQLYRHPGRSQPEFL DAIDRREDTFYVVSFRRDHLPAISHNKTSRPMKSLVMPAMAPNETVSGRGP PGDYEEMMQIECEVMDTRVIHIKTSTVPPSLRKQSPSPGNTTGGPLPGSAAS PAHQASQPLYLNHP
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Atf6b
<b>Protein Names</b>	Recommended name: Cyclic AMP-dependent transcription factor ATF-6 beta Short name= cAMP-dependent transcription factor ATF-6 beta Alternative name(s): Activating transcription factor 6 beta Short name= ATF6-beta cAMP response element-binding p
<b>Expression Region</b>	2-699
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.