



# Recombinant Human Delta-like protein 1 (DLL1)

<b>Product Code</b>	CSB-CF006947HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O00548
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	QVWSSGVFELKLQEFVNKKGLLGNRNCCRGGAGPPPCACRTFFRVLKHYQ ASVSPEPPCTYGS AVTPVLGVDSFSLPDGGGADSAFSNPIRFPFGFTWPGTFS LIIEALHTDSPDDLATENPERLISRLATQRHLTVGEEWSQDLHSSGRTDLKYSY RFVCDEHYGEGCSVFCRPRDDAFGHFTCGERGEKVCNPGWKGPYCTEPIC LPGCDEQHGFCDKPGECKCRVGVWQGRYCDICIRYPGCLHGTCQQPWQCNC QEGWGGLFCNQDLNYCTHHKPKNGATCTNTGQGSYTCSCRPGYTGATCEL GIDECDPSPCKNGGSCDLENSYSCTCPPGFYKICELSAMTCADGPCFNNG RCSDSPDGGYSCRCPVGYSGFNCEKKIDYCSSSPCSNGAKCVDLGDAYLCR CQAGFSGRHCDNVDDCASSPCANGGTCDRGVNDFSCTCPPGYTGRNCSA PVSRCEHAPCHNGATCHERGHRYVCECARGYGGPNCQFLPELPPGPAVVD LTEKLEGQGGPPFWAVCAGVILVLMMLLGC AAVVVCVRLRLQKHRPPADPC RGETETMNNLANCQREKDISVSIIGATQIKNTNKKADFDHSDKNGFKARY PAVDYNLVQDLKGDDTAVRDAHSKRDTKCQPQGSSEEGKTPPTTLRGGEAS ERKRPDSGCSTSKDTKYQSVYVISEEKDECVIATEV
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	DLL1
<b>Protein Names</b>	Recommended name: Delta-like protein 1 Alternative name(s): Drosophila Delta homolog 1 Short name= Delta1 Short name= H-Delta-1
<b>Expression Region</b>	18-723
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.