



Recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 19 (Adam19)

Product Code	CSB-CF001279MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35674
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	EDLHSMKYVELYLVDYAEFQKNRHDQDATKRKLMEIANYVDKFYRSLNIRIAL VGLVWTHGDKCEVSENPYSTLWSFLSWRRKLLAQKSHDNAQLITGRSFQGT TIGLAPLMAMCSVYQSGGVSMDSHSENAIGVASTVAHEIGHNFGMSHDSAHC SASAADGGCIMAATGHPFPKVFWSWNRKELDRYLQTGGGMCLSNMPDTRT LYGGRRCGNGYLEDGEECDCEEEEECKNPCCNASNCTLKEGAECAHGSCC HQCKLVAPGTQCREQVRQCDLPEFCTGKSPHCPTNYYQMDGTPCEGGQAY CYNGMCLTYQECCQQLWGPGARPALDLCFERVNAAGDTYGNCGKGLNGQY RKCSPRDAKCGKIQCCSTQARPLESNAVSIDTTITLNGRRIHCRGTHVYRGPE EEEGEGDMLDPGLVMTGTCKGHNHICFEGQCRNTSFFETEGCGKCKNGHGV CNNKNCHCFPGWSPFCNTPGDGGSVDSGPLPPKSVGPVIAGVFSALFVLA VLVLLCHCYRQSHKLGKPSALPFKLRHQFSCPFVRSQSGGTGHANPTFKLQT PQGKRKVTNTPESLRKPSHPPRPPPDYLRVESPPAPLPAHLNRAAGSSPEA GARIERKESARRPPSRPMPPAPNCLLSQDFSRPRPPQKALPANPVPQGRTG PRSGGTSLLQPPTSGPQPPRPPAVPVPKLPEYRSQRVGAIISSKI
Source	in vitro E.coli expression system
Target Names	Adam19
Protein Names	Recommended name: Disintegrin and metalloproteinase domain-containing protein 19 Short name= ADAM 19 EC= 3.4.24.- Alternative name(s): Meltrin-beta
Expression Region	205-920
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.