



Recombinant Mouse Platelet glycoprotein Ib alpha chain (Gp1ba)

Product Code	CSB-CF009685MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35930
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>QHTCSISKVTSLLEVNCENKKLTALPADLPADTGILHLGENQLGTFSTASLVHFT HLTYLYLDRCELTSLQTNGKLIKLENLDLSHNNLKSLSPLGWALPALTTLDVSN KLGSLSPGVLGSLQLQELYLQNNDLKSLPPGLLLPTTKLKKLNLANNKLRELP SGLLDGLEDLDTLYLQRNWLRTIPKGGFTLLLPFVFLHANSWYCDCEILYFRH WLQENANNVYLWKQGVVDVKDTPNVASVRCANLDNAPVYSYPGKGCPTSSG DTDYDDYDDIPDVPATRTEVKFSTNTKVHTTHWSLLAAPSTSQDSQMISLPP THKPTKKQSTFIHTQSPGFTTLPETMESNPTFYSLKLNTVLIPSPTTLEPTSTQA TPEPNIQPMLTTSTLTTPEHSTTPVPTTTLTTPPEHSTIPVPTTALTTPKPSTIPV PTTATLTTLEPSTTPVPTTATLTTPEPSTTLVPTTATLTTPEHSTTPVPTTATLTT PEHSTTPVPTTATLTTPEPSTTLTNLVSTISPVLTTTTLTTPESTPIETILEQFFTTE LTLPTLESTTTIPEQNSFLNLPEVALVSSDTSSESPFLNSDFCCFLPLGFYVLG LLWLLFASVVLILLTWTWHVTPHSLDMEQSAALATSTHTTSLEVQRARQVTM PRAWLLFLQGSLPTFRSSLFLWVRPNRVRGVLVAGRRPSALSQGRGQDLLGT VGIRYSGHSL</p>
Source	in vitro E.coli expression system
Target Names	Gp1ba
Protein Names	Recommended name: Platelet glycoprotein Ib alpha chain Short name= GP-Ib alpha Short name= GPIb-alpha Short name= GPIbA Short name= Glycoprotein Ibalpha Alternative name(s): CD_antigen= CD42b
Expression Region	17-734
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.