



# Recombinant Human Pro-neuregulin-3, membrane-bound isoform (NRG3)

<b>Product Code</b>	CSB-CF016079HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P56975
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	MSEGAAAASPPGAASAAAASAEEGTAAAAAAAAAGGGPDGGGEGAAEPPRE LRCSDCIVWNRQQTWLCVVPLFIGFIGLGLSLMLLKWIVVGSVKEYVPTDLVDS KGMGQDPFFLSKPSSFPKAMETTTTTTTSTTSPATPSAGGAASSRTPNRISTR TTITRAPTRFPGHRVPIRASPRSTTARNTAAPATVPSTTAPFFSSSTLGSRPPV PGTPTQAMPSPWPTAAYATSSYLHDSTPSWTLSPFQDAASSSSSSSSSATTT TPETSTSPKFHTTTTYSTERSEHFKPCRDKDLAYCLNDGECFVIETLTGSHKHC RCKEGYQGVRCDQFLPKTDSILSDPTDHLGIEFMESEEVYQRQVLSISCIIFGIVI VGMFCAAFYFKSKKQAKQIQEQLKVPQNGKSYSLKASSTMAKSENLVKSHVQ LQNYSKVERHPVTALEKMMESSFVGPQSFPEVPSDRGSQSVKHHRSLSLSSC CSPGQRSGMLHRNAFRRTPPSPRSRLGGIVGPAYQQLEESRIPDQDTIPCQGI EVRKTISHLPIQLWCVERPLDLKYSSSGLKTQRNTSINMQLPSRETNPYFNSLE QKDLVGYSSTRASSVPIIPSVGLEETCLQMPGISEVKS IKWCKNSYSADVNVNS IPVSDCLIAEQQEVKILLETVQE QIRILTDARRSEDYELASVETEDSASENTAFLP LSPTAKSEREAQFVLRNEIQRDSALTK
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	NRG3
<b>Protein Names</b>	Recommended name: Pro-neuregulin-3, membrane-bound isoform Short name= Pro-NRG3 Cleaved into the following chain: 1. Neuregulin-3 Short name= NRG-3
<b>Expression Region</b>	1-720
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.