



Recombinant Human Cadherin-3 (CDH3)

Product Code	CSB-CF005052HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P22223
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	DWVVAPISVPENGGKPPQRLNQLKSNKDRDTKIFYSITGPGADSPPEGVFAV EKETGWLLLNKPLDREEIAKYELFGHAVSENGASVEDPMNISIIVTDQNDHKPK FTQDTFRGSVLEGLVPGTSVMQVTATDEDDAIYTYNGVVAYSIHSQEPKDPHD LMFTIHRSTGTISVISSGLDREKVPEYTLTIQATDMDGDGSTTTAVAVVEILDAN DNAPMFDPPQKYEAHVPENAVGHEVQRLTVTDLDAPNSPAWRATYLIMGGDD GDHFTITTHPESNQGILTRKGLDFEAKNQHTLYVEVTNEAPFVLKLPSTATIV VHVEDVNEAPVFPVPSKVVEVQEGIPTGEPVCVYTAEDPKENQKISYRILRDP AGWLAMPDPSGQVTAVGTLDREDEQFVRNNIYEMVLAMDNGSPPTTGTGT LLLTLIDVNDHGPVPEPRQITICNQSPVRQVLNITDKDLSPHTSPFQAQLTDDSD IYWTAEVNEEGDTVVLVSLKFLKQDQTYDVHLSLSDHGNKEQLTVIRATVCDCH GHVETCPGPWKGGFILPVLGAVLALLFLLLVLVLRKKRRIKEPLLLPEDDTRD NVFYYGEEGGGEEDQDYDITQLHRGLEARPEVLRNDVAPTIPTPMYRPRPA NPDEIGNFIIENLKAANTDPTAPPYDTLLVFDYEGSGSDAASLSSLTSSASDQD QDYDYLNEWGSRFKKLADMYGGGEDD
Source	in vitro E.coli expression system
Target Names	CDH3
Protein Names	Recommended name: Cadherin-3 Alternative name(s): Placental cadherin Short name= P-cadherin
Expression Region	108-829
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.