



Recombinant Mouse Cadherin-3 (Cdh3)

Product Code	CSB-CF005052MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10287
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	EWVMPPIFVPENGGKGFPPQRLNQLKSNKDRGTKIFYSITGPGADSPPEGVFTIE KESGWLLLHMPLDREKIVKYELYGHAVSENGASVEEPMNISIIVTDQNDNPKPF TQDTRGSLVLEGVMPGTSVMQVATDEDDAVNTYNGVVAYSIIHSQEPKEPHD LMFTIHKSTGTISVISSGLDREKVPEYRLTVQATDMDGEGSTTTAEAVVQILDA NDNAPEFEPQKYEAWVPENEVGHEVQRLTVTDLDVPSNPAWRATYHIVGGD DGDHFTITTHPETNQGVLTTKKGLDFEAQDQHTLYVEVTNEAPFAVKLPTATAT VVVHVKDVNEAPVFPVPPSKVIEAQEGISIGELVCIYTAQDPDKEDQKISYTISR PANWLAVDPSGQITAAGILDREDEQFVKNNVYEMVLATDSGNPPTTGTGTL LLTLTDINDHGPIPEPRQIIICNQSPVPQVLNITDKDLSPNSSPFQAQLTHDS DIY WMAEVSEKGD TVALSLKKFLKQD TYDLHLSLSDHGNREQLT MIRATVCD CHG QVFND CPRPWKGGFILPILGAVLALLLLLALLLVRKKR KVKEPLLLPEDD TRD NVFYYGEEGGGEEDQDYDITQLHRGLEARPEVVL RNDVVPTFIPTMYR PRPA NPDEIGNFIIENLKAANTDPTAPPYD SLLVFDYEGSGSDAASLSSLT TSASDQD QDYNYLNEWGSRFKKLADMYGGGEDD
Source	in vitro E.coli expression system
Target Names	Cdh3
Protein Names	Recommended name: Cadherin-3 Alternative name(s): Placental cadherin Short name= P-cadherin
Expression Region	100-822
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.