



Recombinant Rabbit Integrin beta-8 (ITGB8)

Product Code	CSB-CF011892RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26013
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	ENSRCASSHAVSCSECLALGPDCGWCVHEDFISGGPRSERCDIVSNLISKGCP VDSIEYPSVHVTIPSENEVNTQVTPGEVSIQLRPGAAANFMLKIHPLKKYPVDLY YLVDVSASMHNIEKLN SVGNLDRKMAFFSRDFRLGFGSYVDKTVSPYISHP ERIHNCSDYNLDCMPPHGYIHVLSLTENITEFERAVHRQKISGNIDTPEGGFD AMLQAAVCESHIGWRKEAKRLLLVMTDQTSHLALDSKLAGIVVPNDGNCHLRN NVYVKSTTMEHPSLGQLSEKLIDNNINVIFAVQKGQFHWDLLPLLPGTIAGEI ESKAANLNNLVVEAYQKLISEVKVQVESKVPGVYFNVTAICPDGARKLGMEGC SNVTSSDEVLFNVTVTMEKCSVTGGKNYAIKPIGFNETSKIHIHQNCGCECEA SRGGAACKAAEEAPLDSTCPQCQESQCHQEEAQSPSQGCKAHEDQPVCGR GVCVCGKCLCHKMKLGKVYGYKCEKDDFSCPYHHGSLCAGHGECEAGRCQ CFSGWEGDRCQCPAAAQHCVNSKGQVCSGRGTCVCGRCECSDPRSIGRF CEHCPTCPTACSENWNCVQCLPHNLSQAILDQCRTSCASMEQPYVEQASE CFSSPSYLRIFFIIVTFLIGLLKILIRQVILQWNSSKIKSSSDYRVSASKDKLIL QSVCTRAVTYRREKPEEIKLDISKLN AHETFRCNF
Source	in vitro E.coli expression system
Target Names	ITGB8
Protein Names	Recommended name: Integrin beta-8
Expression Region	43-768
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.