



Recombinant Rat P-selectin (Selp)

Product Code	CSB-CF020978RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P98106
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>WTYNYSTKAYSWNNSRAFCKRHFTDLVAIQNKNEIAHLNDVIPYVNSYYWIGIR KINNKWTWVGNTKLTAEAEWADNEPNNKRNNQDCVEIYIKSNSAPGKWND EPCFKRKRALCYTASCQDMSQNSQGERIETIGSYTCSCYPGFYGPCEYVQE CGKFDIPQHVLNCSHPLGDFSFSSQCTFSCPEGYDLNGPSEMQLASGIWT NNPPQCKAVQCQSLEAPLHGTMDCTHPLAAFAYDSSCKFECQPGYRMRGSD ILHCTDSGQWSEPLPTCEAIACEPLESPLHGSMDCFPSTGAFGYNSSCTFRCT EGFVLMGNDAIHCADLGQWTAPAPVCEALQCQEFVPSKAQVSCSDPFGPLK YQSACSFSCDEGSLLVGASVIRCLATGHWSEAPPECQAVSCTPLLSPENGTM TCIQPLGHSNYKSTCQFMCEGDFYLSGPERLDCSPSGHWTGSPPMCEAIKCP EIFAPEQGSLLDCSHVHGEFSVGSTCHFSCNEEFELLGSRNVECTVSGRWSAP PPTCKGVTSLPVPSVRCPALTPGQGTMSCRHHLESFGPNTTCYFGCKTGFT LRGANSRLRCGASGQWTAVTPVCRAVKCSELHMDTAVAMNCSNPWGNFSYG STCAFHCPEGQSLNGSARTTCGEDGHWS DAMPTCQAGTLTIQEALTYLGGAL ASTSGLAVGGTLLALLRKRRLRKKDDGKCP LNPHSHLGTYG VFTNAAYDPTP</p>
Source	in vitro E.coli expression system
Target Names	Selp
Protein Names	Recommended name: P-selectin Alternative name(s): CD62 antigen-like family member P Granule membrane protein 140 Short name= GMP-140 Leukocyte-endothelial cell adhesion molecule 3 Short name= LECAM3 Platelet activation dependent granule-ex
Expression Region	42-768
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.