



Recombinant Chicken Cadherin-1 (CDH1)

Product Code	CSB-CF005034CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08641
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	DWVIPPISCLNHRGPYPMRLVQIKSNKDKEKVVYYSITGQGADSPPVGIFIER ETGWLEVTEQLDREKIDRYTLLSHAVSASGQPVEDPMEIIITVMDQNDNKPVFI KEVFVGYIEENAKPGTSMVTVNATDADDAVNTDNGIVSYSIVSQPPRPHPQM FTIDPAKGIIISVLGTGLDRETPNYTLIVQATDQEGKGLSNTATAIIEVTDANDNI PIFNPTMYEGVVEENKPGTEVARLTVTDQDAPGSPAQQAVYHIKSGNLDGAF SIITDPSTNNGILKTAKGLDYETKSRYDLVVTVENKVPLSVPITLSTASVLVTVLD VNEPPVFVPPIKRVGVPEDLPVGQQVTSYTAQDPDRDMRQKITYRMGSDPAG WLYIHPENGIVTATQPLDRESVHAINSTYKAILAVDNGIPDTTGTGTLTLLLLQDV NDNGPTPEPRSFECISRQPEKQILSIVDKDLPPTYPFKAALEHGSSNNWTVEI RGQDELAMGLKKELEPGEYNIFVKLTDSQGKAQVTQVKAQVCECEGTAKNCE RRSYIVGGLGVPAILGILGGILALLILLILLLLFARRRKVEKEPLLPEDDMRDNV YNYDEEGGGGEEDQDYDLSQLHRGLDARPEVIRNDVAPPLMAAPQYRPRPAN PDEIGNFIDENLKAADTDPTAPPYDSLLVFDYEGGGSEATSLSSLNSSASDQD QDYDYLNEWGNRFKKLAELYGGGEDDE
Source	in vitro E.coli expression system
Target Names	CDH1
Protein Names	Recommended name: Cadherin-1 Alternative name(s): Epithelial cadherin Short name= E-cadherin Liver cell adhesion molecule Short name= L-CAM
Expression Region	161-887
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.