



Recombinant Mouse Stromal interaction molecule 2 (Stim2)

Product Code	CSB-CF022830MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P83093
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	CDLVPRHLRGRRASGSAGAAASPSAAAAGERQALLTDPCMSLSPPCFTEEDR FSLEALQTIHKQMDDDDKGGIEVDESDEFIREDMKYKDANKHSHLHREDKHIT VEDLWKQWKTSEVHNWTLEDTLQWLIEFVELPQYEKNFRDNNVKGTTLPRIA VHETSFMISQLKISDRSHRQKLQKALDVVLFGLTRPPHNWMKDFILTISIVIG VGGCWFAYTQNKTSKEHVAKMMKDLESLQTAEQSLMDLQERLEKAQEENRT VAVEKQNLERKMMDEINYAKEEACRLRELREGAECESRRQYAEQELEQVRM ALKKAEKEFELRSSWSVPDALQKWLQLTHEVEVQYYNIKRQNAEMQLAIKDE AEKIKKKRSTVFGTLHVAHSSSLDEVDHKILEAKKALSELTTCLRERLFRWQQIE KICGFQIAHNSGLPSLTSSLYSDHSWVMPRVSIPPYPIAGGVDDLDDETPPIV PQFPGTVAKPAGSLARSSSLCRSRRSIVPSSPQSQAQLPAHAPLAAHPRHP HHPQHPQHSLPSPDPDILSVSSCPALYRNEEEEEAIYFTAQWEVPDTASEC DSLNSSSSGRKPSPSSLEMYQTLSSRKISRDELSLEDSSRGESPVADSVSRGS PECVGLTETKSMIFSPASRVYNGILEKSCSMHQLSSGIPVPHPRHTSCSSAGN DSKPVQEASNVS RVSSIPHDLCHNGEKSKKPSKIKSLFKKKSK
Source	in vitro E.coli expression system
Target Names	Stim2
Protein Names	Recommended name: Stromal interaction molecule 2
Expression Region	15-746
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.