



Recombinant Mouse Cadherin-9 (Cdh9)

Product Code	CSB-CF005058MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70407
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	GWMWNQFFLLEEYTGTDQYVVGKLHTDQDKGDGNLKYILTGDGAGNLFVIDE NTGDIHAAKRLDREEKSLYLRAKAIDRKTGRQVEPESEFIIKIHDINDNEPKFTK DLYTASVPEMSGVGTSVIQVTATDADDANYGNSAKVVYSILQGQPYFSVDPESE GIIKTALPDMRENKEQYQVVIQAKDMGGQMGGLSGTTTNNITLTDVNNNPPR FPQSTYQFNSLESAPLGLGRIKANDPDMGENAELEYSIAEGEGSDMFDVIT DKDTQEGIIIVKQNLDFEKKMLYTLRVDASNTHPDRFLHLGPFKDSAMVKISV EDVDEPPVFSKLSYLMEVDEVDKEGSIIGQVTAYDPDAMNNIIKYSVDRHTDM DRVFSIHSENGSIFTLKPLDRESSPWHNITITATEINNPKQSSQIPVIFIRILDINDH APEFATYYETFVCENAKSGQLIQTISVMDKDDPPRGHKFFFEFVPEFPLNPFT IVDNKDNTAGIVTRKDGYSRNMNTYLLPVLIFDNDYPIQSSTGTLTIRVCACDN LGNMQSCNAEALMLAAGLSTGALIAILLCVVILLTLIVLFAALKRQRKKEPLISKD DVRDNIVTYNDEGGGEEDTQAFDIGTLRNPEAREDSKLRRDVMPIETIFQIRRT VPLWENIDVQDFIHRRLKENDSDPSAPPYDSLATYAYEGNDSVANSLSLES TADCNQDYDYLSDWGRPFKLAEMYGGNDS
Source	in vitro E.coli expression system
Target Names	Cdh9
Protein Names	Recommended name: Cadherin-9 Alternative name(s): T1-cadherin
Expression Region	53-786
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.