



Recombinant Mouse Cadherin-10 (Cdh10)

Product Code	CSB-CF005035MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70408
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	GWMWNQFFLLEEYTGSDYQYVVGKLSHQDQKGDGSLKYILSGDGAGTLFIIDE KTGDIHATTRIDREEKAFYTLRAQAINRRTL RPVEPESEFVIKIHNDINEPTFPE EIYTASVPEMSVVGTSVVQVTATDADDPSYGN SARVIYSILQGQPYFSVEPET GIIRTALPNMNRENKEQYQVVIQAKDMGGQMGGLSGTTTTVNITLTDVNDNPPR FPQNTIHLRVLESSPVGTAVGSVKATDADTGKNAEVDYRIIDGDGTD MFDIITEK DTQEGIITVKKPLDYENRRLYTLKVEAENTHVDPRFY YLGPFKDTTIVKISIEDVD EPPVFSRSSYLFEVHEDIEVGTIIGTVMARDPDSTSSPIRFTLDRHTDLDRIFNIH SGNGSLYTSKPLDRELSQWHNLT VIAAEINNP KETTRVSVFVRILDVNDNAPQF AVFYDTFVCENARPGQLIQTISAVDKDDPLGGQKFFFSLAAVNPNFTVQDNED NTARILTRKNGFN RHEISTYLLPVVISDNDYPIQSSTGTLTIRVCACDSQGNMQS CSAEALLLPAGLSTGALIAILLCIILLVIVVLF AALKRQRKKEPLILSKEDIRDNI VS YNDEGGGEEDTQAFDIGTLRNPA AIEEKKLRRDIIPETLFI PRRTPTAPDNTDVR DFINERLKEHDL DPTAPPYDSLATYAYEGNDSVAESLSSLES GTTEGDQNYDY LREWGPRFNKLAEMYGGGESDKDA
Source	in vitro E.coli expression system
Target Names	Cdh10
Protein Names	Recommended name: Cadherin-10 Alternative name(s): T2-cadherin
Expression Region	55-788
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.