



Recombinant Rat Cadherin-6 (Cdh6)

Product Code	CSB-CF005055RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55280
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<pre>SWMWNQFFLLEEYTGSDYQYVVGKLSHQDRGDGSLKYILSGDGAGDLFIINE NTGDIQATKRLDREEKPVYILRAQAINRRTGRPEPESEFIIKIHNDNEPIFTK DVYTATVPEMADVGTFFVQVTATDADDPTYGNSAKVVYSILQGQPYFSVESET GIIKTALLNMDRENREQYQVVIQAKDMGGQMGGLSGTTTNNITLTDVNDNPPR FPQSTYQFKTPESPPTPIGRIKASDADVGENAEIEYSITDGEGHDMFDVITD QETQEGIITVKKLLDFEKKRVYTLKVEASNPHEPRFLYLGPFKDSATVRIVVDD VDEPPVFSKPAYILQIREDAQINTTIGSVAQDPDAARNPVKYSVDRHTDMDRI FNIDSGNGSIFTSKLLDRETLWHNITVIATEINNPQSSRVPLYIKVLDVNDNAP EFAEFYETFVCEKAKADQLIQLHAVDKDDPYSGHQFSFLAPEAASGSNFTIQ DNKDNTAGILTRKNGYNRHEMSTYLLPVVISDNDYPVQSSTGTVTVRVCACDH HGNMQSCHAEALIHPTGLSTGALVAILLCIVILLVTVVLFAALRRQRKKEPLIISKE DIRDNIVSYNDEGGGEEDTQAFDIGTLRNPKPWRQQSRRDMVPEALFLPRRT PTARDNTDVRDFISQRLRKMNTDPTAPPYDSLATYAYEGTGSVADSLSSLESV TTDGDQDYGYLSDWGPFRFKKLADMYGGMDSKDS</pre>
Source	in vitro E.coli expression system
Target Names	Cdh6
Protein Names	Recommended name: Cadherin-6 Alternative name(s): Kidney cadherin Short name= K-cadherin
Expression Region	54-789
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.