



Recombinant Mouse Cadherin-6 (Cdh6)

Product Code	CSB-CF005055MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97326
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	SWMWNQFFLLEEYTGSDYQYVVGKLSHQDRGDGSLKYILSGDGAGDLFIINE NTGDIQATKRLDREEKPVYILRAQAVNRRTGRPVEPESEFIIKIHDINDNEPIFTK DVYTATVPEMADVGTFFVQVTATDADDPTYGNSAKVVYSILQGQPYFSVESET GIIKTALLNMDRENREQYQVVIQAKDMGGQMGGLSGTTTNNITLTDVNDNPPR FPQSTYQFKTPESPPTPIGRIKASDADVGENAEIEYSITDGEGHEMFDVITD QETQEGIITVKKLLDFEKKKVYTLKVEASNPVHVEPRFLYLGPFKDSATVRIVVDD VDEPPVFSKLAYILQIREDARINTTIGSVAAQDPDAARNPVKYSVDRHTDMDRIF NIDSGNGSIFTSKLLDRETLWHNITVIATEINNPKQSSRVPLYIKVLDVNDNAPE FAEFYETFVCEKAKADQLIQLRAVDKDDPYSGHQFSFSLAPEAASSSNFTIQD NKDNTAGILTRKNGYNRHEMSTYLLPVVISDNDYPVQSSTGTVTVRVCACDHH GNMQSCHAEALIHPTGLSTGALVAILLCIVILLVTVVLFAALRRQRKKEPLIISKED IRDNIVSYNDEGGGEEDTQAFDIGTLRNPEAMEDSKSRRDIVPEALFLPRRTPT ARDNTDVRDFINQRLKENDTDPTAPPYDSLATYAYEGTGSVADSLSSLESVTT DGDQDYDYLSDWGPFRFKKLADMYGGMDSKDS
Source	in vitro E.coli expression system
Target Names	Cdh6
Protein Names	Recommended name: Cadherin-6 Alternative name(s): Kidney cadherin Short name= K-cadherin
Expression Region	54-790
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.