



Recombinant Human Cadherin-8 (CDH8)

Product Code	CSB-CF005057HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55286
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>GWVWNQMFVLEEFSGPEPILVGR LH TD LDPGSKKIKYILSGDGAGTIFQINDVT GDIHAIKRLDREEKAEYTLTAQAVDWETSKPLEPPSEFIKQVQDINDNAPEFLNG PYHATVPEMSILGTSVTNVTATDADDPVYGNSAKLVYSILEGQPYFSIEPETAIL KTALPNMDREAKEYLVVIQAKDMGGHSGGLSGTTTTLTVTLTDVNDNPPKFAQ SLYHFSVPEDVVLGTAIGRVKANDQDIGENAQSSYDIIDGDGTALFEITSDAQA QDGIIRLRKPLDFETKKSYYTLKVEAANVHIDPRFSGRGPFKDTATVKIVVEDADE PPVFSSPTYLLEVHENAALNSVIGQVTARDPDITSSPIRFSIDRHTDLERQFNIN ADDGKITLATPLDRELSVWHNITIIATEIRNHSQISRVPVAIKVLDVNDNAPEFAS EYEAFLCENGKPGQVIQTVSAMDKDDPKNGHYFLYLLPEMVNPNFTIKKNE DNSLSILAKHNGFNQRQKQEVYLLPIIISDSGNPPLSSTSTLTIRVCGCSNDGVVQ SCNVEAYVLPGLSMGALIAILACIILLVIVVLFVTLRRHKNEPLIHKDDEDVRENII RYDDEGGGEEDTEAFDIATLQNP D G I N G F L P R K D I K P D L Q F M P R Q G L A P V P N G VDVDEFINVRLHEADNDPTAPPYDSIQIYGYEGRGVSAGSLSSLESTTS D S D Q NFDYLSDWGPRFKRLGELYSVGESDKET</p>
Source	in vitro E.coli expression system
Target Names	CDH8
Protein Names	Recommended name: Cadherin-8
Expression Region	62-799
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>