



Recombinant Mouse Cadherin-11 (Cdh11)

Product Code	CSB-CF005036MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55288
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	GWVWNQFFVIEEYTGPDPLVGRGLHSDIDSGDGNIKYILSGEGAGTIFVIDDKS GNIHATKTLTREERAQYTLMAQAVDRDNRPLEPPSEFIVKVQDINDNPPEFLH EIIHANVPERSNVGTSVIQVTASDADDPTYGNSAKLVYSILEGQPYFSVEAQTG IIRTALPNMDREAKEEYHVVIQAKDMGGHMGGLSGTTKVTITLTDVNDNPPKFP QSVYQMSVSEAAVPGEEVGRVKAKDPDIGENGLVTYNIVDGDGIELFEITTDYE TQDGVVKKLPVDFETKRAYSLKIEAANVHIDPKFISNGPFKDTVTVKISVEDAD EPPMFLAPSYIHEVQENAAAGTVVGRVHAKDPDAANSPIRYSIDRHTDLDRFFT INPEDGFIKTTKPLDREETAWLNISVFAAEIHNRHQETKVPVAIRVLDVNDNAPK FAAPYEGFICESDHPKALSNOPIVTVSADDQDDTANGPRFIFSLPPEIMHNPNF TVRDNRDNTAGVYARRGGFSRQKQDFYLLPIVISDGGIPPMSSNTLTIKVCGC DVNGALLSCNAEAYILNAGLSTGALIAILACIVILLVIVVLFVTLRRQKKEPLIVFEE EDVRENIITYDDEGGGEEDTEAFDIATLQNPDGINGFIPRKDIKPEYQYMPRPG LRPAPNSVDVDDFINTRIQEADNDPTAPPYDSIQIYGYEGRGSVAGSLSSLESA TTDSLDYDYLQNWGPRFKKLADLYGSKDTFDDDS
Source	in vitro E.coli expression system
Target Names	Cdh11
Protein Names	Recommended name: Cadherin-11 Alternative name(s): OSF-4 Osteoblast cadherin Short name= OB-cadherin
Expression Region	54-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.