



Recombinant Mouse Cadherin-like protein 26 (Cdh26)

Product Code	CSB-CF005050MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59862
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	HQPLHRSKRRWVLTTLLELQEEDPGPFKLVGELFNNMSNNVSLIYLIRGPGVD EFPEIGLFSIEDHQSGKIYVHRPVDREVTPSFMVHFDVDRSTGKVVDESLIFNI RIRDVNDHAPQFPEKEFNISVKESQAAGQPIFQLLTVDLQDENTPNSQVLYFLV SQTPLLRESGFRIDLISGEVRLSGCLHYETAPLFTLIVRASDCGEPSTLSTATHV SVEDSNNHMPTFMEDRYEIRISEGQVEQGVLYLPVQDGDSPFTPAWRTQFNI WDGNEEGHFDIVTDPETNQGLLSIIKPLDYESQVAHSLVVVENQEQLFVCTE GQLQPLRKAMASTMVSQVLDTNDPPAFHPQSFIVSEEDGAWPAIQLGYFNA TDPDRADSQIRYKLVHDPENWVTVDKSGVVTTKKQIDRESPHVNSFYTIIVH AVDNGLPPLTGTGTLMLFLSDVNDNAPTLRPHSRHLEVCEASGSQPLIEAED ADLDPYADPFTFDLNAQGDVEETWMLRTKQEGHSAELTMLRSVPPGDYLV PLFIADRQGLAQKQTVHVRICSCRSGSECEEPSDTWLLWWALSPVGAALMVL SAALLCLLRCSCTFGPKRLRGFIPSDSGHQTLIVYNEESKVPASAQGCDDTFEPR GVKTLSSSTPVYLDRMVPRQQPLQLLEGRVVEAWSQKLQSIDVLEGDTGYLP HVYREEGECEGAETLSSLTFLEQDLSPKLLGCSGSKSTPSEAMCFTSRVPS
Source	in vitro E.coli expression system
Target Names	Cdh26
Protein Names	Recommended name: Cadherin-like protein 26
Expression Region	21-763
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.