



Recombinant Human Glutamate carboxypeptidase 2 (FOLH1)

Product Code	CSB-CF008782HU(A4)
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q04609
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MWNLLHETDSAVATARRPRWLCAGALVLAGGGFFLLGFLFGWFIKSSNEATNIT PKHNMKAFLDELKAENIKKFLYNFTQIPHLAGTEQNFQLAKQIQSQWKEFGLDS VELAHYDVLLSYPNKTHPNYISIINEDGNEIFNTSLFEPPPPGYENVSDIVPPFSA FSPQGMPEGDLVYVNYARTEDFFKLERDMKINCSGKIVIARYGKVFRGNKVKN AQLAGAKGVILYSDPADYFAPGVKSYPDGWNLPGGGVQQRGNILNLNGAGDPL TPGYPANEYAYRRGIAEAVGLPSIPVHPIGYYYDAQKLEKMGGSAPPDSSWRG SLKVPYNVGPFTGNFSTQKVKMHIHSTNEVTRIYNVIGTLRGAVEPDRYVILG GHRDSWVFGGIDPQSGAAVVHEIVRSFGTLKKEGWRPRRTILFASWDAEEFG LLGSTEWAEENSRLQERGVAYINADSSIEGNYTLRVDCTPLMYSLVHNLTKEL KSPDEGFEGKSLYESWTKKSPSPEFSGMPRISKLGSGNDFEVFFQRLGIASGR ARYTKNWETNKFSGYPLYHSVYETYELVEKFYDPMFKYHLTVAQVRGGMVFE LANSIVLPFDCRDYAVVLRKYADKIYSISMKHPQEMKTVSFSDFSLFSAVKNFT EIASKFSERLQDFDKSNPIVLRMMNDQLMFLERAFIDPLGLPDRPFYRHVIYAP SSHNKYAGESFPGIYDALFDIESKVDPSKAWGEVQRQIYVAAFTVQAAAETLSE VA
Source	in vitro E.coli expression system
Target Names	FOLH1
Protein Names	Recommended name: Glutamate carboxypeptidase 2 EC= 3.4.17.21 Alternative name(s): Cell growth-inhibiting gene 27 protein Folate hydrolase 1 Folylpoly- gamma-glutamate carboxypeptidase Short name= FGCP Glutamate carboxypeptidase II Short n
Expression Region	1-750
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.