



Recombinant Human Cadherin-15 (CDH15)

Product Code	CSB-CF005040HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55291
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	LPYPLVQIKSDKQQLGSVIYSIQGPGVDEEPRGVFSIDKFTGKVFLNAMLDREK TDRFRLRAFALDLGGSTLEDPTDLEIVVVDQNDNRPAFLQEAFTRVLEGAVP GTYYVTRAEATDADDPETDNAALRFSILQQGSPELFSIDELTGEIRTVQVGLDRE VVAVYNLTLQVADMSGDGLTATASAIITLDDINDNAPEFTRDEFFMEAIEAVSGV DVGRLEVEDRDLPGSPNWWARFTILEGDPDGQFTIRTDPKTNEGVLISIVKALD YESCEHYELKVSQNEAPLQAAALRAERGQAKVRVHVQDTNEPPVFQENPLR TSLAEGAPPGTLVATFSARDPDTEQLQRLSYSKYDPEDWLQVDAATGRIQT QHVLSASPFLKGGWYRAIVLAQDDASQPRATGTLSIEILEVNDHAPVLAPP PGSLCSEPHQGPGLLLGATDEDLPPHGAPFFHQLSPRLPELGRNWSLSQVNV SHARLRPRHQVPEGLHRLSLLLRDSGQPPQREQPLNVTVCRCGKDGVCPLP GAAALLAGGTGLSLGALVIVLASALLLVLLVALRARFWKQSRGKGLLHGPQ DDL RDNLNYDEQGGGEEEDQDAYDISQLRHPTALSPLGPPPLRRDAPQGRL HPQPPRVLPTSPLDIADFINDGLEAADS DPSVPPYDTALIYDYEGDGSVAGTLS SILSSQGDEDQDYDYL RDWGPFRFARLADMYGHPCGLEYGARWDHQAREGLS PGALLPRHRGRTA
Source	in vitro E.coli expression system
Target Names	CDH15
Protein Names	Recommended name: Cadherin-15 Alternative name(s): Cadherin-14 Muscle cadherin Short name= M-cadherin
Expression Region	61-814
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.