



Recombinant Human Glycosyltransferase-like protein LARGE1 (LARGE)

| | |
|--------------------------|---|
| Product Code | CSB-CF012749HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O95461 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sequence | MLGICRGRRKFLAASLSLLCIPAITWIYLFSGSFEDGKPVSLSPLESQAHSPTYT ASSQRERESLEVRMREVEEENRALRRQLSLAQGRAPSHRRGNHSKTYSMEE GTGDSENLRAGIVAGNSSECGQQPVVEK CETIHVAIVCAGYNASRDVTLVKS VLFHRRNPLHFHFIADSI AEQILATLFQ TWMVPAVRVDFYNAD ELKSEVSWIPN KHYSGIYGLMKLVLT KTLPANLERVIVLDTDITFATDIAELWAVFHKFKGQQV LG LVENQSDWYLG NLWKNHRPWPALGRGYNTGVILLLLDKLRKMKWEQMWRLT AERELMGMLSTSLADQDIFNAVIKQNPFLVYQLPCFWNVQLSDHTRSEQCYR DVSDLKVIHWNSPKKLRVKNKHVEFFRNLYLTFLEYDGNLLRRELFGCPSEAD VNSENLQKQLSELDEDDLCYEFRRERFTVHRTHLYFLHYEYEPAADSTDVTLV AQLSMDRLQMLEAICKHWE GPISLALYLSDAEAQQFLRYAQGSEVLMSRHNV GYHIVYKEGQFYPVNLLRNVAMKHISTPYMFLSDIDFLPMYGLYEYLRKSVIQL DLANTKKAMIVPAFETLRYRLSFPKSKAELLSMLDMGTLFTFRYHVWTKGHAP TNFAKWRTATTPYRVEWEADFE PYVVRRDCPEYDRRFVGF GWNKVAHIME LDVQEYEFIVLPNAYMIHMPHAPSFDITKFRSNKQYRICLKT LKEEFQQDMSRR YGFAALKYLTAENNS |
| Source | in vitro E.coli expression system |
| Target Names | LARGE1 |
| Protein Names | Recommended name: Glycosyltransferase-like protein LARGE1 EC= 2.4.-.- Alternative name(s): Acetylglucosaminyltransferase-like 1A |
| Expression Region | 1-756 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |