



Recombinant Human Integrin beta-3 (ITGB3)

Product Code	CSB-CF011885HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P05106
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	GPNICTTRGVSSCQQCLAVSPMCAWCSDEALPLGSPRCDLKENLLKDNCAPESIEFPVSEARVLEDRPLSDKSGDSSQVTQVSPQRIALRLRPDDSKNFSIQVRQVEDYPVDIYYLMDLSYSMKDDLWSIQNLGTKLATQMRKLTSLNRIGFGAFVDKPVSPYMYISPPEALENPCYDMKTTCLPMFGYKHVLTLDQVTRFNEEVKKQSVSRNRDAPEGGFDAIMQATVCDEKIGWRNDASHLLVFTTDAKTHIALDGRLAGIVQPNDGQCHVGSDNHYSASTTMDYPSLGLMTEKLSQKNINLIFAVTENVVNL YQNYSELIPGTTVGVLSMDSSNVLQLIVDAYGKIRSKVELEVRDLPEELSLSFNATCLNNEVIPGLKSCMGLKIGDTVFSFSIEAKVRGCPQEKEKSFTIKPVGFKDSLIVQVTFDCDCACQAQAEPNSHRCNNGNGTFCGVCRCGPGWLGSGQCECSEEDYRPSQQDECSPPREGQPVCSSQRGECLCGQCVCHSSDFGKITGKYCECDDFS CVRYKGEMCSGHGQCSCGDCLCDSWTGYCNCCTTRTDTCMSSNGLLCSGRGKCECGSCVCIQPGSYGDTCEKCPTCPDACTFKKECVECKKFDRGALHDEN TCNRYCRDEIESVKELKDTGKDAVNCTYKNEDDCVVRVFQYYEDSSGKSILYVV EEPECPKGPDILVVLLSVMGAILLIGLAALLIWKLLITIHDRKEFAKFEEERARAK WDTANNPLYKEATSTFTNITYRGT
Source	in vitro E.coli expression system
Target Names	ITGB3
Protein Names	Recommended name: Integrin beta-3 Alternative name(s): Platelet membrane glycoprotein IIIa Short name= GPIIIa CD_antigen= CD61
Expression Region	27-788
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.