



# Recombinant Human High affinity nerve growth factor receptor (NTRK1)

<b>Product Code</b>	CSB-CF016133HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04629
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	AAPCPDACCPHGSSGLRCTR DGALDSLHHLPGAENL TELYIENQQHLQHLELR DLRGLGELRNLTIVKSGLRFVAPDAFHFTPRLSRLNLSFNALESLSWKT VQGLS LQELVLSGNPLHCSCALRWLQRWEEEG LGGVPEQKLQCHGQGPLAHMPNAS CGVPTLKVQVPNASVDVGDDVLLRCQVEGRGLEQAGWILTELEQSATVMKSG GLPSLGLTLANVTSDLNRKNVTCWAENDVGRAEVSVQVNVSF PASVQLHTAV EMHHWCIPFSVDGQPAPSLRWLFNGSVLNETSFIFTEFLEPAANETVRHGCLR LNQPETHVNNGN YTLAANPFGQASASIMAAFMDNPFEPEDPIPVSFSPVDT NSTSGDPVEKKDETPFGVSVAVGLAVFACLFLSTLLLVLNKCGR RNKFGINRP AVLAPEDGLAMSLHFMTLGGSSLSPT EGKGSGLQGHIIENPQYFSDACVHHIK RRDIVLKWELGEGAFGKVFLAECHNLLPEQDKMLVAVKALKEASE SARQDFQ REALLTMLQHQHIVRFFGVCTEGRPLL MVFEYMRHGD LNRFLRSHGPDALL AGGEDVAPGPLGLGQLLAVASQVAAGMVYLAGLHFVHRDLATRNCLVGGQLV VKIGDFGMSRDIYSTDYRVGGRTMLPIR WMPPE SILYRKFTTESDVWSFGVV LWEIFTY GKQPWYQLSNTEAIDCITQGRELERPRACPPEVYAIMRGCWQREP QQRHSIKDVHARLQALAQAPPVYLDVLG
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	NTRK1
<b>Protein Names</b>	Recommended name: High affinity nerve growth factor receptor EC= 2.7.10.1 Alternative name(s): Neurotrophic tyrosine kinase receptor type 1 TRK1- transforming tyrosine kinase protein Tropomyosin-related kinase A Tyrosine kinase receptor Tyro
<b>Expression Region</b>	33-796
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.