



Recombinant Human Endothelin-converting enzyme-like 1 (ECE1)

Product Code	CSB-CF007373HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95672
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEPPYSLTAHYDEFQEVKYVSRGAGGARGASLPPGFPLGAARSATGARSGL PRWNRREVCLLSGLVFAAGLCAILAAMLALKYLGPVAAGGGACPEGCPERKA FARAARFLAANLDASIDPCQDFYSFACGGWLRRAIPDDKLTGTIAAIGEQNE ERLRLLARPGGGPGGAAQRKVRAFFRSCCLDMREIERLGPRPMLEVEDCGG WDLGGAEERPGVAARWDLNRLLYKAQGVYSAALFSLTVSLDDRNSSRYVIRI DQDGLTLPERTLYLAQDEDESEKILAAAYRVFMERVLSLLGADAVEQKAQEILQVE QQLANITVSEHDDLRRDVSSMYNKVTLGQLQKITPHLRWKWLLDQIFQEDFSE EEEVLLATDYMQQVSQLIRSTPHRVLHNYLVWRVVLSEHLSPPFREALHE LAQEMEGSDKPQELARVCLGQANRHFGMALGALFVHEHFAASKAKVQQLV EDIKYLGRLEELDWMDAETRAAARAKLQYMMVMVGYPDFLLKPDVAVDKEY EFEVHEKTYFKNILNSIRFSIQLSVKKIRQEVKSTWLLPPQALNAYYLPNKNQM VFPAGILQPTLYDPDFPQSLNYGGIGTIIGHELTHGYDDWGGQYDRSGNLLHW WTEASYSRFLRKAECIVRLYDNFTVYNQRVNGKHTLGENIADMGGKLAYHAY QKWWREHGPEHPLPRLKYTHDQLFFIAFAQNWCIKRRSQQSIYLQVLTDKHAPE HYRVLGVSVSQFEFGRGFHCPKDSMPMNAHKCSVW
Source	in vitro E.coli expression system
Target Names	ECEL1
Protein Names	Recommended name: Endothelin-converting enzyme-like 1 EC= 3.4.24.- Alternative name(s): Xce protein
Expression Region	1-775
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.