



Recombinant Human Fibroblast growth factor receptor 4 (FGFR4)

Product Code	CSB-CF008648HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P22455
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>LEAEEVELEPCLAPSLEQQEQELTVALGQPVRLCCGRAERGGHWYKEGSRL APAGRVRGWRGRLEIASFLPEDAGRYLCLARGSMIVLQNLTITGDSLTSND DEDPKSHRDP SNRHSYPQQAPYWTHPQRMEKKLHAVPAGNTVKFRCPAAGN PTPTIRWLKDGQAFHGENRIGGIRLRHQHWSLVMESVVPDRGTYTCLVENA VGSIRYNYLLDVLERSPHRPILQAGLPANTTAVVGSDEVLLCKVYSDAQPHIQW LKHIVINGSSFGADGFPYVQVLKTADINSSEVEVLYLRNVS AEDAGEYTCLAGN SIGLSYQSAWLTVLPEEDPTWTA APEARYTDIILYASGSLALAVLLLLAGLYRG QALHGRHPRPPATVQKLSRFPLARQFSLESGSSGKSSSSLVRGVRLSSSGPA LLAGLVSLDPLDPLWEFPRDRLVLGKPLGEGCFGQVRAEAFGMDPARPDQ ASTVAVKMLKDNASDKDLADLVSEMEVMKLI GRHKNIINLLGVCTQEGPLYVIV ECAAKGNLREFLRARRPPGPDLS PDGPRSSEGPLSFPVLVSCAYQVARGMQY LESRKCIHRDLAARNVLVTEDNVMKIADFGLARGVHHIDYKKT SNGRLPVKW MAPEALFDRVYTHQSDVWSFGILLWEIFTLGGSPYPGIPVEELF SLLREGHRM DRPPHCPELYGLMRECWHAAPSQRPTFKQLVEALDKVLLAVSEEYLDLRLTF GPYSPSGGDASSTCSSSDSVFSDPLPLGSSSFPFGSGVQT</p>
Source	in vitro E.coli expression system
Target Names	FGFR4
Protein Names	Recommended name: Fibroblast growth factor receptor 4 Short name= FGFR-4 EC= 2.7.10.1 Alternative name(s): CD_antigen= CD334
Expression Region	22-802
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.