



Recombinant Mouse Protein FAM171A2 (Fam171a2)

Product Code	CSB-CF008131MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2A699
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	KSPPEPPSPQEILIKVQVYVSGELVPLARASVDVFGNRTLLAAGTTDSEGVATL PLSYRLGTWVLVTAARPGFLTNSVPWRVDKLPLYASVSLYLLPERPATLILYED LVHILLGSPGARSQPWVQFQRRARLPVSSTYSQLWASLTPASTQQEMRAFP AFLGTEASSSGNGSWLELIPLAAVSVHLLTGNNGTEVPLSGPIHLSLPVSEPR LAVGTSIPAWRFDPKSGLWVRNGTG VIRKEGRQLYWTFVSPQLGYWVAAMA SPTSGLVTITSGIQDIGTYHTIFLLTILAAALALLVLILLCLLIYYCRRRCLKPRQQHR KLQLSGPSDNKRQATSMSQLHLICGGPLEPTSSGDPEAPPPGSLHSAFSSS RDLASSRDDFFRAKPRSASRPAAEPPGARTVEGAGLKSARSVEGPGGLEPSL DEYRRGPAGAAFLHEPPSPPPSFDHYLGHKGAAESKTPDFLLSQSVDQLAR PPSLSQPGQLIFCGSIDHLKDNVYRNVMP TLVIPAHYVRLGGEAGAAGVGDEA TPPEGSAAGPARPFPQDPQRPLMQGHAGAGGDSGGGEGWGGGGRSAPVS GSVTIPVLFNESTMAQLNGELQALTEKKLLELGVKPHPRAWFVSLDGRSNSQV RHSYIDLQAGSGGRSTDASLD SGVDVHEARPARRRLPREERERAQLPAPPPP APPRLALSEDTEPSSSESRTGLCSPEDNSLTPLLDEVVAPEGRAATVPRGRGR SRGDSSRSSASELRRDSLTSPEDELGAEVGDEAGDKKSPWQRREERPLMVF NVK
Source	in vitro E.coli expression system
Target Names	Fam171a2
Protein Names	Recommended name: Protein FAM171A2
Expression Region	30-822
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.