



Recombinant Human Protein FAM171A2 (FAM171A2)

Product Code	CSB-CF008131HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8MVW0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	KSPPEPPSPQEILIKVQVYVSGELVPLARASVDVFGNRTLLAAGTTDSEGVATL PLSYRLGTWVLVTAARPGFLTNSVPWRVDKLPLYASVSLYLLPERPATLILYED LVHILLGSPGARSQPLVQFQRRARLPVSSTYSQLWASLTPASTQQEMRAFPA FLGTEASSSGNGSWLELMPLTAVSVHLLTNGTEVPLSGPIHLSLPVSETRAL TVGTSIPAWRFDPKSGLWVRNGTG VIRKEGRQLYWTFVSPQLGYWVAAMAS PTAGLVTITSGIQDIGTYHTIFLLTILAAALLLVILLCLLIYYCRRRCLKPRQQHRK LQLSGPSDGNKRDQATSMSQLHLICGGPLEPAPSGDPEAPPPGPLHSAFSSS RDLASSRDDFFRTKPRSASRPAAEPSGARGGESAGLKGARSAEGPGGLEPG LEEHRGSPSGAAAFLEPPSPPPFDHYLGHKGAEGKTPDFLLSQSVDQLA RPPSLGQAGQLIFCGSIDHLKDNVYRNVMP TLVIPAHYVRLGGEAGAAGVGDE PAPPEGTAPGPARAFPQDPQRPQMPGHSGPGGEGGGGGEGWGAGRAA PVSGSVTIPVLFNESTMAQLNGELQALTEKKLLELGVKPHPRAWFVSLDGRSN SQVRHSYIDLQAGGGARSTDASLD SGVDVHEARPARRRPAREERERAPPAAP PPPPAPPRLALSEDTEPSSSESRTGLCSPEDNSLTPLLDEVAAPEGRAATVPR GRGRSRGDSSRSSASELRRDSLTSPEDELGAEVGDEAGDKKSPWQRREERP LMVFNVK
Source	in vitro E.coli expression system
Target Names	FAM171A2
Protein Names	Recommended name: Protein FAM171A2
Expression Region	30-826
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.