



Recombinant Chicken Fibroblast growth factor receptor 2 (FGFR2)

Product Code	CSB-CF008645CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P18461
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	RPSFNLVVEDATLEPEEPPTYQISQPDVHSALPGEPELELRCLKDAVMISWT KDGVPPLGPDNRTVIIGEYLQIKDASPRDSGLYACTAIRTLDSDTLYFIVNVTDAL SSGDEDDNDGSEDFVNDNSQMRAPYWHTDKMEKRLHAVPAANTVKFRCP AMGNPTPTMRWLKNGKEFKQEHRIGGYKVRNQHWLIMESVVPSPDKGNYTCI VENQYGSINHTYHLDVVERSHPRPILQAGLPANASAVVGGDVEFVCKVYSDAQ PHIQWIKHVERNGSKYGPDGLPYLQVLKAAGVNTTDEIEVLYIRNVTFEDAGE YTCLAGNSIGISFHTAWLTVLPAPEKEKEFPTSPDYLEIAIYCIGVFLIACMVLTVI LCRMKNTTKKPDFSSQPAVHKLTKRIPLRRQVTVSADSSSSMNSNTPLVRITT RLSSTADAPMLAGVSEYELPEDPKWEFPRDKLTLGKPLGEGCFGQVMAEAV GIDKDRPKEAVTVAVKMLKDDATEKDLSDLVSEMEMMKMIGKHKNIINLLGACT QDGPLYVIVEYASKGNLREYLRARRPPGMEYSFDINRVPEEQMTFKDLVSCY QLARGMEYLASQKCIHRDLAARNVLVTENNVMKIADFLARDINNIDYKKTNN GRLPVKWMPEALFDRVYTHQSDVWSFGVLMWEIFTLGGSPYPGIPVEELFK LLKEGHRMDKPANCTNELYMMMRDCWQAVPSQRPTFKQLVEDLDRILTLTNN EEYLDLSGPLEQYSPSYDTRSSCSSGDDSVFSPDPMPYEPCLPKYQHMNGS VKT
Source	in vitro E.coli expression system
Target Names	FGFR2
Protein Names	Recommended name: Fibroblast growth factor receptor 2 Short name= FGFR-2 EC= 2.7.10.1 Alternative name(s): Tyrosine kinase receptor CEK3
Expression Region	24-823
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.