



Recombinant Human Fibroblast growth factor receptor 2 (FGFR2)

Product Code	CSB-CF008645HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21802
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	RPSFSLVEDTTLEPEEPPTKYQISQPEVYVAAPGESLEVRCLLKDAAVISWTKD GVHLGPNNRTVLIGEYLQIKGATPRDSGLYACTASRTVDSETWYFMVNVDAI SSGDEDDTDGAEDFVSENSNNKRAPYWTNTEKMEKRLHAVPAANTVKFRC PAGGNPMPMRWLKNGKEFKQEHRIGGYKVRNQHWSLIMESVVPDCKGNYT CVVENEYGSINHTYHLDVVERSHPRPILQAGLPANASTVVGGDVEFVCKVYSD AQPHIQWIKHVEKNGSKYGPDGLPYLKVLAAGVNTTDKEIEVLYIRNVTFEDA GEYTCLAGNSIGISFHSAWLTVLPAPGREKEITASPDYLEIAIYCIGVFLIACMVV TVILCRMKNNTTKPDFSSQPAVHKLTKRIPLRRQVTVSAESSSSMNSNTPLVRI TTRLSSTADTPMLAGVSEYELPEDPKWEFPRDKLTLGKPLGEGCFGQVMAE AVGIDKDKPKEAVTVAVKMLKDDATEKDLSDLVSEMEMMKMIGKHKNIINLLGA CTQDGPLYVIVEYASKGNLREYLRRARRPPGMEYSYDINRVPEEQMTFKDLVSC TYQLARGMEYLASQKCIHRDLAARNVLVTENNVMKIADFGLARDINNIDYYKKT TNGRLPVKWMapeALFDRVYTHQSDVWSFGVLMWEIFTLGGSPYPGIPVEEL FKLLKEGHRMDKPANCTNELYMMMRDCWHAVPSQRPTFKQLVEDLDRILT TNEEYLDLSQPLEQYSPSYDPDRSSCSSGDDSVFSPDPMPYEPCLPQYPHIN GSVKT
Source	in vitro E.coli expression system
Target Names	FGFR2
Protein Names	Recommended name: Fibroblast growth factor receptor 2 Short name= FGFR-2 EC= 2.7.10.1 Alternative name(s): K-sam Short name= KGFR Keratinocyte growth factor receptor CD_antigen= CD332
Expression Region	22-821
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.