



Recombinant Mouse CD63 antigen (Cd63)

Product Code	CSB-CF004950MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41731
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	AVEGGMKCVKFLLYVLLLAFCACAVGLIAIGVAVQVVLKQAITHETTAGSLLPVV IIAVGAFRLFVAFVGGCGACKENYCLMITFAIFLSLIMLVEVAVAIAGYVFRDQVK SEFNKSFQQMQNYLKDNKTATILDKLQKENNCCGASNYTDWENIPGMAKDR VPDSCCINITVGCNDFKESTIHTQGCVETIAIWLKRNILLVAAAALGIAFVEVLGI IFSCCLVKSIRSGYEVV
Source	in vitro E.coli expression system
Target Names	Cd63
Protein Names	Recommended name: CD63 antigen Alternative name(s): CD_antigen= CD63
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.