



Recombinant Guinea pig B2 bradykinin receptor (BDKRB2)

Product Code	CSB-CF002652GU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70526
Product Type	Transmembrane Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Sequence	MFNITSQVSALNATLAQGN SCLDAEWW SWLNTIQAPFLWVLFVLAVLENIFVLS VFFLHKSSCTVAEIYLG NLAVADLILAFGLPFWAITIANNFDWLFGEVLCRMVNT MIQMNMYS SIFLMLVSI DRYLALVKTMSMGRMRGVRWAKLYSLVIWGCALLL SSPMLVFRTMKDYRDEGHNVTA CLIIYPSLTWQVFTNVLLNLVGFLLPLSIITFC TVQIMQVLRNNEMQKFKEIQTERRATV LVLAVLLL FVVCWLPFQIGTFLDTRLRL GFLPGCWEHVIDLITQISSYLAYSNSCLNPLVYVIVGKRFRKKSREVYHGLCRS GGCVSEPAQSENSMGT LR T S I S V D R Q I H K L Q D W A R S S S E G T P P G L L
Source	in vitro E.coli expression system
Target Names	BDKRB2
Protein Names	Recommended name: B2 bradykinin receptor Short name= B2R Short name= BK-2 receptor
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.