



Recombinant Chicken Frizzled-7 (FZD7)

Product Code	CSB-CF009110CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O57329
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	AHHEDKAISVPDHGFCQPISIP LCTDIAYNQ TILPNLLGHTNQEDAGLEVHQFYPLVKVQCSAELKFFLCSMYAPVCTVLEQAIPPCRSLCERARQGCEALMNKFGFQWPERLRCENFPVHGAGEICV GQNTSDAPP GPGGAGGRGATAQPTAGYLPDLTPPQPAAGFSFSCPRQLKVPPYLG YRFLGERDCGAPCEPGRPNGLMYFKEAEVRFARLWVGWVSVLCCASTLFTVLT YLVDMRRFSYPERPIIFLSGCYFMVAVAYAAGFLL EERVVCLERFSEDGYRTVAQGTKKEGCTILF MILYFFGMASSIWWVILSLTWFLAAGMKWGHEAIEANSQYFH LAAWAVPAVKTITILAMGQVDGDVLSGVCYVGIYSVDSL RGFVLAPLFVYLFIGTSFLLAGFVSLFR IRTIMKHDGKTEKLEKLMVRIGVFSVLYTVPATIVVAC YFYEQAFRSTWEKTWLLQTCKTYAVPCP SHFAPMSPDFTVFM IKYLMTMIVGITTGF WIWWSGKTLQSWRR FYHRLSTGSKGETAV
Source	in vitro E.coli expression system
Target Names	FZD7
Protein Names	Recommended name: Frizzled-7 Short name= Fz-7 Short name= cFz-7
Expression Region	32-567
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.