



Recombinant Mouse Annexin A4 (Anxa4), partial

Product Code	CSB-EP001845MO
Relevance	Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.
Abbreviation	Recombinant Mouse Anxa4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P97429
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	EAKGGTVKAASGFNATEDAQTLRKAMKGLGTDEDIIGILAYRNTAQRQEIRSA YKSTIGRDLIEDLKSELSSNFEQVILGLMTPTVLYDVQELRRAMKGAGTDEGCLI EILASRTPEEIRRINQTYQQQYGRSLEEDICSDFMFQRVLSLSAAGRDEGN YLDDALMKQDAQELYEAGEKRWGTDEVKFLSILCSRNRNHLHVFDEYKRISQ KDIEQSIKSETSGSFEDALLAIVKCMRSKPSYFAERLYKSMKGLGTDDNTLIRV MVSRAEIDMLDIRASFKRLYGKSLYSFIKGDTSGDYRKVLLVLCGGDD
Research Area	Cardiovascular
Source	E.coli
Target Names	Anxa4
Protein Names	Annexin IV Annexin-4 Synonyms:Anx4
Expression Region	2-319aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	40.8 kDa
Protein Length	Partial
Image	

