




# Recombinant Human Protein Dr1 (DR1)

<b>Product Code</b>	CSB-EP007175HU
<b>Relevance</b>	The association of the DR1/DRAP1 heterodimer with TBP results in a functional repression of both activated and basal transcription of class II genes. This interaction precludes the formation of a transcription-competent complex by inhibiting the association of TFIIA and/or TFIIB with TBP. Can bind to DNA on its own. Component of the ATAC complex, a complex with histone acetyltransferase activity on histones H3 and H4.
<b>Abbreviation</b>	Recombinant Human DR1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q01658
<b>Alias</b>	Down-regulator of transcription 1 Negative cofactor 2-beta
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	ASSSGNDDDLTIPRAAINKMIKETLPNVRVANDARELVVNCCTEFIHLLISSEANEI CNKSEKKTISPEHVIALESLSLGFSGYISEVKEVLQECKTVALKRRKASSRLENL GIPEEELLRQQQELFAKARQQQAELAQQEWLQMQQAAQQAQLAAASASASN QAGSSQDEEDDDDI
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Source</b>	E.coli
<b>Target Names</b>	DR1
<b>Protein Names</b>	Recommended name: Protein Dr1 Alternative name(s): Down-regulator of transcription 1 Negative cofactor 2-beta Short name= NC2-beta TATA-binding protein-associated phosphoprotein
<b>Expression Region</b>	1-176aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	46.3kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.