




Recombinant Mouse Ornithine transcarbamylase, mitochondrial (Otc)

Product Code	CSB-EP017270MO
Abbreviation	Recombinant Mouse Otc protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11725
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	SQVQLKGRDLLTLKNFTGEEIQYMLWLSADLKFRKQKGEYLPLLQGKSLGMIF EKRSTRTRLSTETGFALLGGHPSFLTQQDIHLGVNESLTDARVLSSMTDAVLA RVYKQSDLDLTAKEASIPVINGLSPLYHPIQILADYLTQEHYGSLLKGLTSLWIG DGNNILHSIMMSAAKFGMHLQAATPKGYEPDPNIVKLAEQYAKENGTKLSMTN DPLEAARGGNVLITDTWISMGQEDEKKKRLQAFQGYQVTMKTAKVAASDWF LHCLPRKPEEVDDEVFYSRSLVFPPEAENRKWTIMAVMVSLLDYSPVLQKPK F
Research Area	Others
Source	E.coli
Target Names	Otc
Expression Region	33-354aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	52.1kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.