




# Recombinant Human Thyroid hormone receptor alpha (THRA)

<b>Product Code</b>	CSB-EP023510HU
<b>Relevance</b>	Isoform Alpha-1: Nuclear hormone receptor that can act as a repressor or activator of transcription. High affinity receptor for thyroid hormones, including triiodothyronine and thyroxine.
<b>Abbreviation</b>	Recombinant Human THRA protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P10827
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MEQKPSKVECGSDPEENSARSPDGKRKRKNGQCSLKTSMSGYIPSYLDKDE QCVVCGDKATGYHYRCITCEGCKGFFRRTIQKNLHPTYSCKYDSCCVIDKTR NQCQLCRFKKCIAVGMAMDLVLDDSKRVAKRKLIEQNRERRRKEEMIRSLQQ RPEPTPEEWDLIHATEAHRSTNAQGSHWKQRRKFLPDDIGQSPIVSMFDGDK VDLEAFSEFTKIITPAITRVVDFAKKLPMFSELPCEQIILLKGCCMEIMSLRAAV RYDPESDTLTLSGEMAVKREQLKNGGLGVVSDAIFELGKLSAFNLDDTEVAL LQAVLLMSTDRSGLLCVDKIEKSQEAYLLAFEHYVNRKHNIPHFWPKLLMKE REVQSSILYKAAAEGRPGGSLGVHPEGQQLLGMHVVQGPQVRQLEQQQLGE AGSLQGPVLQHQSPKSPQQRLELLHRSGILHARAVCGEDDSSEADSPSSSE EEPEVCEDLAGNAASP
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	THRA
<b>Protein Names</b>	Nuclear receptor subfamily 1 group A member 1 V-erbA-related protein 7 Short name: EAR-7 c-erbA-1 c-erbA-alpha EAR7, ERBA1, NR1A1, THRA1, THRA2
<b>Expression Region</b>	1-490aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-B2M-tagged
<b>Mol. Weight</b>	68.8 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .