




Recombinant Human Transcobalamin-1 (TCN1)

Product Code	CSB-EP023317HU
Relevance	Binds vitamin B12 with femtomolar affinity and protects it from the acidic environment of the stomach.
Abbreviation	Recombinant Human TCN1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20061
Alias	Haptocorrin Short name: HC Protein R Transcobalamin I Short name: TC I Short name: TCI
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	EICEVSEENYIRLKPLLNTMIQSNYNRGTSAVNVVLSLKLVGIIQITLMQMIIQQI KYNVKSRLSDVSSGELALIILALGVCRNAEENLIYDYHLIDKLENKFQAEIENME AHNGTPLTNYYQLSLDLALCLFNGNYSTAEEVNHFTPENKNYYFGSQFSVDT GAMAVLALTCVKKSLINGQIKADEGSLKNISIYTKSLVEKILSEKKENGLIGNTF TGEAMQALFVSSDYNNENDWNCQQTTLNLTVEISQGAFSNPNAQAQVLPALM GKTFLDINKDSSCVSASGNFNISADEPITVTPPDSQSYISVNYSVRINETYFTNV TVLNGSVFLSVMEKAQKMNDTIFGFTMEERSWGPYITCIQGLCANNNDRTYW ELLSGGEPLSQGAGSYVVRNGENLEVRWSKY
Research Area	Signal Transduction
Source	E.coli
Target Names	TCN1
Protein Names	Recommended name: Transcobalamin-1 Short name= TC-1 Alternative name(s): Transcobalamin I Short name= TC I Short name= TCI
Expression Region	24-433aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	61.6kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.