




Recombinant Human B-lymphocyte antigen CD19 (CD19), partial

Product Code	CSB-EP004888HU
Relevance	Assembles with the antigen receptor of B-lymphocytes in order to decrease the threshold for antigen receptor-dependent stimulation.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15391
Alias	B-lymphocyte surface antigen B4 Differentiation antigen CD19 T-cell surface antigen Leu-12
Product Type	Recombinant Protein
Immunogen Species	Recombinant Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	PEEPLVVKVEEGDNAVLQCLKGTSDGPTQQLTWSRESPLKPFLKLSLGLPGLG IHMRPLAIWLFIFNVSQQMGGFYLCQPGPPSEKAWQPGWTVNVEGSGELFR WNVSDLGGLGCGLKNRSSEGPSSPSGKLMSPKLYVWAKDRPEIWEGEPPCL PPRDSL NQSL SQDLTMAPGSTLWLSCGVPPDSVSRGPLSWTHVHPKGP KSL L SLELKDDRPARDMWVMETGLLLPRATAQDAGKYYCHRG NLTMSFHLEITARP VLWHWLLRTGGWK
Research Area	Immunology
Source	E.coli
Target Names	CD19
Expression Region	20-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	34.1kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.