




Recombinant Mouse Tryptophan 5-hydroxylase 1 (Tph1)

Product Code	CSB-EP024100MO
Relevance	L-tryptophan + tetrahydrobiopterin + O ₂ = 5-hydroxy-L-tryptophan + 4a-hydroxytetrahydrobiopterin.
Abbreviation	Recombinant Mouse Tph1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17532
Alias	Tryptophan 5-monoxygenase 1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MIEDNKENKENDHSSERGRVTLIFSLNEVGGLIKVLKIFQENHVSL LHIESRK SKQRNSEFEIFVDCDISREQLNDIFPLLKSHATVLSVDSPDQLTAKEDVMETVP WFPKKISDLDFCANRVLLYGSELDADHPGFKDNVYRRRRKYFAELAMNYKHG DPIPKIEFTEEEIKTWGTIFRELNKLYPTHACREYLRNLP LLSKYCGYREDNIPQL EDVSNFLKERTGFSIRPVAGYLSRDFLSGLAFRVFHCTQYVRHSSDPLYTPE PDTCHELLGHVPLLAEPSFAQFSQEIGLASLGASEETVQKLATCYFFTVEFGLC KQDQQLRVFGAGLLSSISELKHALSGHAKVKPFDPKIACKQECLITSFQDVYFV SESFEDAKEKMREFAKTVKRPFGLKYNPYTQSVQVLRDTSITSAMNELRYDL DVISDALARVTRWPSV
Research Area	Signal Transduction
Source	E.coli
Target Names	Tph1
Expression Region	1-447aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	67.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.