



# Recombinant Mouse Lysosomal acid Glucosylceramidase (Gba1)

<b>Product Code</b>	CSB-MP009289MO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P17439
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	A QPCIPKSFYGY SSVVCVCNAS YCDSLDPVTL PALGTFSRYE STRRGRRMEL SVGAIQANRT GTGLLLTLQP EKKFQKVKGF GGAMTDATAL NILALSPPTQ KLLRSYFST NGIEYNIIRV PMASCDFSIR VYTYADTPND FQLSNFSLPE EDTKLKIPLI HQALKMSSRP ISLFASPWTS PTWLKTNGRV NGKGSCLKGQP GDIFHQTWAN YFVKFLDAYA KYGLRFWAVT AENEPTAGLF TGYPFQCLGF TPEHQRFDIS RDLGPALANS SHDVKLLMLD DQRLLLPWA EVVLSDEAA KYVHGIAVHW YMDFLAPAKA TLGETHRLFP NTMLFASEAC VGSKFWEQSV RLGSWDRGMQ YSHSIITNLL YHVTGWDWN LALNPEGGPN WVRNFVDSPI IVDIPKDAFY KQPMFYHLGH FSKFIPEGSQ RVALVASEST DLETVALLRP DGSVVVVLN RSSEVDPLTI SDPDLGFLET VSPGYSIHTY LWRRQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Gba1
<b>Protein Names</b>	Recommended name: Glucosylceramidase EC= 3.2.1.45 Alternative name(s): Acid beta-glucosidase Beta-glucocerebrosidase D-glucosyl-N-acylsphingosine glucohydrolase
<b>Expression Region</b>	20-515
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.