




Recombinant Mouse Proto-oncogene Wnt-3 (Wnt3)

Product Code	CSB-EP026135MO
Relevance	Ligand for mbers of the frizzled family of seven transmbane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.
Abbreviation	Recombinant Mouse Wnt3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17553
Alias	Proto-oncogene Int-4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GYPIWWSLALGQQYTSLASQPLLCSIPGLVPKQLRFCRNYIEIMPSVAEGVKL GIQECQHQFRGRRWNCTTIDDSLAIFGPVLDKATRESAFVHAIASAGVAFVTR SCAEGTSTICGCDSHHKGPPGEGWKWGGCSEDADFGVLVSREFADARENRP DARSAMNKHNNNEAGRRTILDHMHKCKCHGLSGSCEVKTCWWAQPDFRAIG DFLKDKYDSASEMVVEKHRESRGWVETLRAKYALFKPPTERDLVYYENSPNF CEPNPETGSFGTRDRTCNVTSHGIDGCDLLCCGRGHNTRTEKRKEKCHCVFH WCCYVSCQECIRIYDVHTCK
Research Area	Others
Source	E.coli
Target Names	Wnt3
Protein Names	Recommended name: Proto-oncogene Wnt-3 Alternative name(s): Proto-oncogene Int-4
Expression Region	22-355aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	53.5kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.