




Recombinant Rabbit Integrin beta-8 (ITGB8), partial

Product Code	CSB-EP011892RB
Relevance	Integrin alpha-V/beta-8 is a receptor for fibronectin.
Abbreviation	Recombinant Rabbit ITGB8 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26013
Product Type	Recombinant Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PVDLYYLVDVSASMHNIEKLNLSVGNLSRKMAFFSRDFRLGFGSYVDKTVSP YISIHPERIHNQCSDYNLDCMPPHGYIHVLSLTENITEFERAVHRQKISGNIDTP EGGFDAMLQAAVCESHIGWRKEAKRLLLVMTDQTSHLALDSKLAGIVVPNDG NCHLRNNVYVKSTTMEHPSLGQLSEKLIDNNINVIFAVQGKQFHWYKDLLPLLP GTIAGEIESKAANLNNLVVEAYQKL
Research Area	Signal Transduction
Source	E.coli
Target Names	ITGB8
Protein Names	Recommended name: Integrin beta-8
Expression Region	146-384aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	31.9kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.