




Recombinant Human Melanoma-associated antigen 1 (MAGEA1)

Product Code	CSB-EP013323HU
Abbreviation	Recombinant Human MAGEA1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43355
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSLEQRSLHCKPEEALAEQAQAEALGLVCVQAATSSSSPLVLGTLEEVPTAGSTD PPQSPQGASAFPTTINFTRQRQPSEGSSSREEEGPSTSCILESIFRAVITKKVA DLVGFLLLKYRAREPVTKAEMLESVIKNYKHCFPEIFGKASESLQLVFGIDVKEA DPTGHSYVLVTCLGLSYDGLLDGNQIMPKTGFLIIVLVMIAAMEGGHAPEEEIWE ELSVMEVYDGREHSAYGEPRKLLTQDLVQEKYLEYRQVPDSDPARYEFLWGP RALAETSYVKVLEYVIKVSARVRFFFPSLREAALREEEEGV
Research Area	Cancer
Source	E.coli
Target Names	MAGEA1
Protein Names	Recommended name: Melanoma-associated antigen 1 Alternative name(s): Antigen MZ2-E Cancer/testis antigen 1.1 Short name= CT1.1 MAGE-1 antigen
Expression Region	1-309aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	40.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.