





Recombinant Mouse Protein Wnt-7b (Wnt7b)


Product Code	CSB-EP026142MO
Relevance	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.
Abbreviation	Recombinant Mouse Wnt7b protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28047
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ALSSVVALGANIICNKIPGLAPRQRAICQSRPDAIIVIGEGAQMGIDECQHQRFRF GRWNC SALGEKTVFGQELRVGSREAAFTYAITAAGVAHAVTAACSQGNLSNC GCDREKQGYYNQAEGWKWGGCSADVRYGIDFSRRFVDAREIKKNARRLMNL HNNEAGRKVLDRMKLECKCHGVSGSCTTKTCWTTLPKFREVGHLLKEKYNA AVQVEVVRASRLRQPTFLRIKQLRSYQKPMETDLVYIEKSPNYCEEDAATGSV GTQGRLCNRTSPGADGCDTMCCGRGYNTHQYTKVWQCNCCKFWCCFVKC NTCSERTEVFTCK
Research Area	Stem Cells
Source	E.coli
Target Names	Wnt7b
Protein Names	Wnt-7b
Expression Region	25-349aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Mol. Weight	56.3 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP026142MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) Wnt7b.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP026142MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) Wnt7b.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.